

ESYSTA[®]

● Digital Diabetes Care

ESYSTA[®] Bluetooth[®] Low Energy Insulin Pen Profile

Version 04.00

Emperra[®]
E-Health Technologies

Emperra[®] GmbH
E-Health Technologies

Zeppelinstraße 48a
D-14471 Potsdam
Germany

Phone: +49 (0) 331 - 979 34 80 - 0

Fax: +49 (0) 331 - 979 34 80 - 19

info@emperra.com
www.emperra.com

Table of contents

- 1 General 4
 - 1.1 Scope 4
 - 1.2 Roles 4
 - 1.3 Roles and Relationships 5
 - 1.4 Requirements 5
- 2 Services and Characteristics 6
 - 2.1 Overview 6
 - 2.2 Insulin Pen Service 6
 - 2.2.1 Insulin Dose Characteristic 6
 - 2.2.2 Insulin Dose Context Characteristic 10
 - 2.2.3 Insulin Pen Feature Characteristic 10
 - 2.2.4 RACP Characteristic 12
 - 2.3 Insulin Pen Custom Service 14
 - 2.3.1 Insulin Pen Custom Value 14

1 General

This specification summarizes the data transfer protocol used for the communication between the ESYSTA[®] BT Pen and the ESYSTA[®] App for mobile devices. The whole data transfer operates over a Bluetooth Low Energy transport only.

The transfer procedure is similar to the official Bluetooth Generic Attribute Profile (GATT)-based Glucose Profile¹ and Glucose Service². For further information regarding used terms, abbreviations, data formats etc. please refer to these resources.

Note: All fields are shown in the order of LSO to MSO, where LSO = Least Significant Octet and MSO = Most Significant Octet.

1.1 Scope

This specification is designed to be suitable as an official standard for transferring data from any Insulin Pen that is equipped with a Bluetooth Low Energy interface. It thus covers features that might be supported by future devices, but that are currently not supported by the ESYSTA BT Pen.

1.2 Roles

This document defines two roles: Insulin Pen and Collector. The Insulin Pen is the device that is used to inject insulin and the Collector is the device that receives the Insulin Pen application data and other related data from an Insulin Pen.


At any given time, an Insulin Pen shall be connected to only one Collector.

¹ Glucose Profile, <https://www.bluetooth.com/specifications/specs/glucose-profile-1-0/>

² Glucose Service, <https://www.bluetooth.com/specifications/specs/glucose-service-1-0/>

1.3 Roles and Relationships

The following diagram shows the relationships between service and profile roles.



Note: Profile roles are represented by blue boxes and the services are represented by grey boxes

An Insulin pen instantiates the Insulin Pen Service and the Device Information Service. The Collector instantiates the Current Time Service.

1.4 Requirements

The Collector shall instantiate the Current Time Service with Local Time Information, as described in the Current Time Profile (CTP).

2 Services and Characteristics

2.1 Overview

Service	UUID	
Generic Access	00001800-0000-1000-8000-00805f9b34fb	
Generic Attribute	00001801-0000-1000-8000-00805f9b34fb	
Device Information	0000180a-0000-1000-8000-00805f9b34fb	
Battery Information	0000180f-0000-1000-8000-00805f9b34fb	
Insulin Pen	c2f00001-1466-4478-b6dc-88c865226d5f	
Insulin Pen Custom	c2f00c01-1466-4478-b6dc-88c865226d5f	

2.2 Insulin Pen Service

To receive insulin dose values, the Collector needs to use the RACP.

All insulin dose entries are delivered through the Insulin Dose characteristic and optional with the Insulin Dose Context characteristics. The Collector needs to subscribe the Client Characteristic Configuration for “Insulin Dose”, “Insulin Dose Context” and “Record Access Control Point” characteristics.

Characteristic	Assigned Number	Properties
Insulin Dose	c2f0000d-1466-4478-b6dc-88c865226d5f	Notify
Insulin Dose Context	c2f0000c-1466-4478-b6dc-88c865226d5f	Notify
Insulin Pen Feature	c2f0000f-1466-4478-b6dc-88c865226d5f	Read
RACP	00002a52-0000-1000-8000-00805f9b34fb	Indicate, Write

2.2.1 Insulin Dose Characteristic

The “Insulin Dose” is used to deliver all base information for a given entry using the following structure. Each entry is delivered through a notification.

Name	Field Requirement	Format	Additional Information																																																												
Flags	Mandatory	8bit	<table border="1"> <thead> <tr> <th>Bit</th> <th>Name</th> <th>Key</th> <th>Value</th> <th>Requires</th> </tr> </thead> <tbody> <tr> <td rowspan="2">0</td> <td>Time Offset Present</td> <td>0</td> <td>False</td> <td></td> </tr> <tr> <td></td> <td>1</td> <td>True</td> <td>C.1</td> </tr> <tr> <td rowspan="2">1</td> <td>Local Time Present</td> <td>0</td> <td>False</td> <td></td> </tr> <tr> <td></td> <td>1</td> <td>True</td> <td>C.2</td> </tr> <tr> <td rowspan="2">2</td> <td rowspan="2">Insulin Dose, Type and Injection Location Present</td> <td>0</td> <td>False</td> <td></td> </tr> <tr> <td>1</td> <td>True</td> <td>C.3</td> </tr> <tr> <td rowspan="2">4</td> <td rowspan="2">Insulin Dose Unit</td> <td>0</td> <td>IU</td> <td></td> </tr> <tr> <td>1</td> <td>L</td> <td></td> </tr> <tr> <td rowspan="2">5</td> <td rowspan="2">Device Status Annunciation Present</td> <td>0</td> <td>False</td> <td></td> </tr> <tr> <td>1</td> <td>True</td> <td>C.4</td> </tr> <tr> <td>6</td> <td colspan="4">reserved for further use</td> </tr> <tr> <td rowspan="2">7</td> <td rowspan="2">Context Information Follows</td> <td>0</td> <td>False</td> <td></td> </tr> <tr> <td>1</td> <td>True</td> <td></td> </tr> </tbody> </table>	Bit	Name	Key	Value	Requires	0	Time Offset Present	0	False			1	True	C.1	1	Local Time Present	0	False			1	True	C.2	2	Insulin Dose, Type and Injection Location Present	0	False		1	True	C.3	4	Insulin Dose Unit	0	IU		1	L		5	Device Status Annunciation Present	0	False		1	True	C.4	6	reserved for further use				7	Context Information Follows	0	False		1	True	
			Bit	Name	Key	Value	Requires																																																								
			0	Time Offset Present	0	False																																																									
					1	True	C.1																																																								
			1	Local Time Present	0	False																																																									
					1	True	C.2																																																								
			2	Insulin Dose, Type and Injection Location Present	0	False																																																									
					1	True	C.3																																																								
			4	Insulin Dose Unit	0	IU																																																									
					1	L																																																									
5	Device Status Annunciation Present	0	False																																																												
		1	True	C.4																																																											
6	reserved for further use																																																														
7	Context Information Follows	0	False																																																												
		1	True																																																												
Sequence Number	Mandatory	uint16																																																													
Base Time	Mandatory																																																														
Year		unit16																																																													
Month		unit8																																																													
Day		unit8																																																													
Hour		unit8																																																													
Minute		unit8																																																													
Second		unit8																																																													
Time Offset	C.1	sint16																																																													
Local Time	C.2																																																														
Time Zone		sint8																																																													
DST Offset		uint8																																																													
Insulin Dose	C.3	sfloat																																																													
Dose Type	C.3	nibble	<table border="1"> <thead> <tr> <th>Key</th> <th>Values</th> </tr> </thead> <tbody> <tr> <td>0</td> <td>Reserved for future use</td> </tr> <tr> <td>1</td> <td>Regular injection</td> </tr> </tbody> </table>	Key	Values	0	Reserved for future use	1	Regular injection																																																						
Key	Values																																																														
0	Reserved for future use																																																														
1	Regular injection																																																														

			2	Priming defined by user																																						
			3	Priming enforced by pen after power up																																						
			4	Priming enforced by pen the carpule was changed																																						
			5	Priming enforced by pen because the cannula was change																																						
			6	Priming enforced by pen because of an unclear device state																																						
			15	Insulin application type not available																																						
Injection Location	C.3	nibble	<table border="1"> <thead> <tr> <th>Key</th> <th>Values</th> </tr> </thead> <tbody> <tr> <td>0</td> <td>Reserved for future use</td> </tr> <tr> <td>1</td> <td>Abdomen</td> </tr> <tr> <td>2</td> <td>Buttocks</td> </tr> <tr> <td>3</td> <td>Hips</td> </tr> <tr> <td>4</td> <td>Upper arm</td> </tr> <tr> <td>5</td> <td>Lower arm</td> </tr> <tr> <td>6</td> <td>Upper leg</td> </tr> <tr> <td>7</td> <td>Lower leg</td> </tr> <tr> <td>15</td> <td>Injection location not available</td> </tr> </tbody> </table>		Key	Values	0	Reserved for future use	1	Abdomen	2	Buttocks	3	Hips	4	Upper arm	5	Lower arm	6	Upper leg	7	Lower leg	15	Injection location not available																		
Key	Values																																									
0	Reserved for future use																																									
1	Abdomen																																									
2	Buttocks																																									
3	Hips																																									
4	Upper arm																																									
5	Lower arm																																									
6	Upper leg																																									
7	Lower leg																																									
15	Injection location not available																																									
Device Status Annunciation	C.4	int16	<table border="1"> <thead> <tr> <th>Bit</th> <th>Name</th> <th>Key</th> <th>Value</th> </tr> </thead> <tbody> <tr> <td rowspan="2">0</td> <td rowspan="2">Device battery low at time</td> <td>0</td> <td>False</td> </tr> <tr> <td>1</td> <td>True</td> </tr> <tr> <td rowspan="2">1</td> <td rowspan="2">Device lifetime exceeded at time</td> <td>0</td> <td>False</td> </tr> <tr> <td>1</td> <td>True</td> </tr> <tr> <td rowspan="2">2</td> <td rowspan="2">Device lifetime near at time</td> <td>0</td> <td>False</td> </tr> <tr> <td>1</td> <td>True</td> </tr> <tr> <td rowspan="2">3</td> <td rowspan="2">Entered dose higher than the device can deliver</td> <td>0</td> <td>False</td> </tr> <tr> <td>1</td> <td>True</td> </tr> <tr> <td rowspan="2">4</td> <td rowspan="2">Entered dose lower than the device can deliver / negative dose</td> <td>0</td> <td>False</td> </tr> <tr> <td>1</td> <td>True</td> </tr> <tr> <td>5</td> <td></td> <td>0</td> <td>False</td> </tr> </tbody> </table>		Bit	Name	Key	Value	0	Device battery low at time	0	False	1	True	1	Device lifetime exceeded at time	0	False	1	True	2	Device lifetime near at time	0	False	1	True	3	Entered dose higher than the device can deliver	0	False	1	True	4	Entered dose lower than the device can deliver / negative dose	0	False	1	True	5		0	False
Bit	Name	Key	Value																																							
0	Device battery low at time	0	False																																							
		1	True																																							
1	Device lifetime exceeded at time	0	False																																							
		1	True																																							
2	Device lifetime near at time	0	False																																							
		1	True																																							
3	Entered dose higher than the device can deliver	0	False																																							
		1	True																																							
4	Entered dose lower than the device can deliver / negative dose	0	False																																							
		1	True																																							
5		0	False																																							

				Application Procedure was interrupted	1	True
			6	Application Procedure timed out	0	False
					1	True
			7	A carpule change preceded the application	0	False
					1	True
			8	The carpule was missing at the time of the application	0	False
					1	True
			9	A carpule of wrong type was inserted at the time of the application	0	False
					1	True
			10	A cannula change preceded the application	0	False
					1	True
			11	The cannula was missing at the time of the application	0	False
					1	True
			12	The device temperature was too high since the last carpule change	0	False
					1	True
			13	The device temperature was too low since the last carpule change	0	False
					1	True
			15	Time fault has occurred and time may be inaccurate	0	False
			14	reserved for further use		

2.2.2 Insulin Dose Context Characteristic

The “Insulin Dose Context” is only used to deliver concrete temperature values for a given time. Each entry is delivered through a notification.

Name	Field Requirement	Format	Additional Information				
Flags	Mandatory	8bit	Bit	Name	Key	Value	Requires
			0	Temperature Value Present	0	False	
					1	True	C.1
			1	Temperature Unit	0	Celsius	
					1	Fahrenheit	
			2	Device Error Event Present	0	False	C.2
1	True						
Sequence Number	Mandatory	uint16					
Temperature Value	C.1	float					
Device Error Event	C.2	uint8	Bit	Name	Key	Value	
			0	General device fault has occurred at time	0	False	
					1	True	
			1	Bluetooth transmission fault has occurred at time	0	False	
					1	True	
			2	Device sensor error at time	0	False	
					1	True	
3-7	reserved for further use						

2.2.3 Insulin Pen Feature Characteristic

The “Insulin Pen Feature” characteristic is used to determine all features supported by device.

Name	Field Requirement	Format	Additional Information				
Flags	Mandatory	8bit	Bit	Name	Key	Value	Requires
			0	Low Battery Detection Supported	0 1	False True	
			1	End of Life Detection Supported	0 1	False True	
			2	Priming Detection Supported	0 1	False True	
			3	Dose High-Low Detection Supported	0 1	False True	
			4	Insulin Temperature High-Low Detection Supported	0 1	False True	
			5	Application Interrupt Detection Supported	0 1	False True	
			6	Application Interrupt Detection Supported	0 1	False True	
			7	Carpule Change Detection Supported	0 1	False True	
			8	Cannula Change Detection Supported	0 1	False True	
			9	Insulin Type Detection Supported	0 1	False True	
			10	Insulin Brand Detection Supported	0 1	False True	

			11	General Device Fault Supported	0	False			
					1	True			
			12	Time Fault Supported	0	False			
					1	True			
			13	Multiple Bond Supported	0	False			
					1	True			
			14	Carpule Error Detection Supported	0	False			
					1	True			
			15	Carpule Error Detection Supported	0	False			
					1	True			

2.2.4 RACP Characteristic

The “Record Access Control Point” characteristic provides basic management functionality for the Insulin Pen record database. This enables functions including counting records, transmitting records and clearing records based on filter criterion.

Name	Field Requirement	Format	Additional Information														
Op Code	Mandatory	uint8	<table border="1"> <thead> <tr> <th>Key</th> <th>Description</th> </tr> </thead> <tbody> <tr> <td>0</td> <td>Reserved for future use (Operator: N/A)</td> </tr> <tr> <td>1</td> <td>Report stored records</td> </tr> <tr> <td>3</td> <td>Abort Operation</td> </tr> <tr> <td>4</td> <td>Report number of stored records</td> </tr> <tr> <td>5</td> <td>Number of stored records response</td> </tr> <tr> <td>6</td> <td>Response Code</td> </tr> </tbody> </table>	Key	Description	0	Reserved for future use (Operator: N/A)	1	Report stored records	3	Abort Operation	4	Report number of stored records	5	Number of stored records response	6	Response Code
			Key	Description													
			0	Reserved for future use (Operator: N/A)													
			1	Report stored records													
			3	Abort Operation													
			4	Report number of stored records													
			5	Number of stored records response													
6	Response Code																
Operator	Mandatory	uint8	<table border="1"> <thead> <tr> <th>Key</th> <th>Description</th> </tr> </thead> <tbody> <tr> <td>0</td> <td>Null</td> </tr> <tr> <td>1</td> <td>All records</td> </tr> <tr> <td>3</td> <td>Greater than or equal to</td> </tr> </tbody> </table>	Key	Description	0	Null	1	All records	3	Greater than or equal to						
			Key	Description													
			0	Null													
			1	All records													
3	Greater than or equal to																

			5	First record (i.e. oldest record).
			6	Last record (i.e. most recent record).
Operand Filter Type	Mandatory	uint8	Key Description	
			0	N/A
			1	Filter Sequence Number
Operand Filter Parameter	Mandatory	variable	Filter operand, e.g. sequence number	
Operand Response Code	Mandatory	uint8	Key Description	
			0	N/A
			1	Success
			2	Op Code not supported
			3	Invalid Operator (received operator does not meet the requirements of service)
			4	Operator not supported
			5	Invalid Operand (received operand does not meet the requirements of service)
			6	No records found
			7	Abort unsuccessful
			8	Procedure not completed
			9	Operand not supported

2.3 Insulin Pen Custom Service

To terminate the connection, the Collector needs to write a 1 to the value.

Characteristic	Assigned Number	Properties
Insulin Dose Custom Value	c2f00c0a-1466-4478-b6dc-88c865226d5f	Read, Write

2.3.1 Insulin Pen Custom Value

Name	Field Requirement	Format	Additional Information
Value	Mandatory	8bit	When a Collector no longer needs the connection, is shall write a value of 1 into this field to terminate the connection and save power on both devices.